Miltenyi Biotec

gentleMACS™ Dissociator

User manual

Version 04
Original instructions

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. MACS is a registered trademark, gentleMACS is a trademark of Miltenyi Biotec GmbH.

Copyright © 2010 Miltenyi Biotec GmbH. All rights reserved.

No part of this user manual may be reproduced, stored in a retrieval system, transmitted, published, or distributed in any form or by any means, electronically, mechanically, by photocopying, microfilming, recording, or otherwise, without the prior written consent of Miltenyi Biotec GmbH; however, notwithstanding the foregoing, the owners of the gentleMACS™ Dissociator may make copies solely for purposes of training personnel in the use and servicing of the unit within their business or organization.

Maximal care has been taken by Miltenyi Biotec GmbH in the preparation of this user manual. However, Miltenyi Biotec GmbH shall not be liable for any technical or editorial errors or omissions contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this document. The information in this document is provided "as is" without warranty of any kind, and is subject to change without notice.

gentleMACS™

Dissociator

User manual

Miltenyi Biotec GmbH
Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany

Phone +49-2204-8306-0
Fax +49-2204-85197
macstec@miltenyibiotec.de

www.miltenyibiotec.com

Version 04
Original instructions

Contents

Chapter 1	7	Important information
	11	Wichtige Informationen
	15	Informations importantes
	19	Informazioni importanti
	23	Información importante
Chapter 2	27	Introduction
Chapter 3	28	gentleMACS™ Products and Spare Parts
Chapter 4	30	Storage and handling
Chapter 5	31	Installation
Chapter 6	32	Operating the gentleMACS™ Dissociator
	32	General information
	33	Description of the instrument
	33	Layout of the keyboard
	33	Panel buttons
	34	gentleMACS™ Program Card
	34	Function of the LED
	34	Display
	34	Power ON
	35	Function of the menus
	35	Select and run a program
	36	Copy a program
	37	Delete a program
	37	Emergency shutdown procedure
Chapter 7	38	Performance and limitations of use
Chapter 8	39	Maintenance
	39	Cleaning
	40	Cleaning procedure
	40	Disassembly for cleaning and reassembly
	41	Safety
	42	Internal quality control
	42	Disposal information
Chapter 9	43	Troubleshooting
Chapter 10	44	Technical specifications
Chapter 11	45	Limited warranty
Chapter 12	46	Technical service
Chapter 13	47	EC Declaration of conformity

Please read before use!

Please read all information contained in this chapter before use. The use of an instrument for a task other than that for which it is intended, as well as improper or insufficient maintenance can greatly reduce the life time of the instrument and will invalidate the instrument's warranty. Incorrect handling and care as well as misuse can lead to premature wear or can cause hazards to users. The following chart is an illustrated glossary depicting the symbols that are used in this user manual and on the gentleMACS™ Dissociator.

Symbols

Warning

Triangles indicate general warnings. Warnings include the risk of damage to the equipment, severe personal injury, or loss of life. Consult the user manual for further instructions and proceed with caution.

Warning

This symbol indicates the risk of loss of life or severe injury to the instrument operator due to high voltage. Consult the user manual for further instructions and proceed with caution.

Warning

This symbol indicates the risk of loss of life, severe injury to the instrument operator, or equipment damage due to potentially dangerous chemical or biological material. Consult the user manual for further instructions and proceed with caution.

Caution

This symbol indicates the risk of equipment damage if instructions in the user manual are not followed precisely. It is therefore strongly recommended that the instrument operator consult the corresponding chapters in the respective documentation.

**Protective conductor terminal**

Symbol is attached on the inside of the instrument. Warning for service personnel only.

**On (supply)****Off (supply)**

Warnings and precautions

This chapter contains information for the safe operation of the gentleMACS™ Dissociator. Detailed instructions regarding safety are provided throughout this manual, marked with the symbols described above.

The operator is responsible for the correct use of the gentleMACS Dissociator following the instructions given in this manual and in full accordance with all local, state, and federal laws that may apply. The impeccable and safe operation of the instrument presumes that any operator understands the user manual and pays attention to all safety information.

Note: All safety instructions in this user manual must be strictly adhered to.

Keep this user manual in a safe place, accessible for anyone using the gentleMACS Dissociator.

At all times, local working area safety instructions, laboratory policies, and standards regarding laboratory safety and health, and prevention of accidents must be adhered to.

In the case of safety-relevant defects of the equipment, Miltenyi Biotec technical service must be contacted.

To ensure safe operation of the gentleMACS Dissociator all covers or parts that have been removed from the gentleMACS Dissociator must be re-installed prior to setting the instrument into operation. Do not allow fluids or objects to enter the interior of the instrument.

Incorrect handling of the gentleMACS Dissociator bears risks in the form of mechanical, electrical, chemical, and biological hazards.

Electrical hazards

Please be particularly aware of electrical hazards. Standard electrical safety precautions must be applied at all times. Only qualified personnel are permitted to perform electrical servicing. It is strictly forbidden to operate the instrument with disassembled or damaged parts.

The instrument is equipped with a three-wire electrical grounding-type plug. This is a safety feature. Do not insert the plug into a non-grounded power outlet. If you cannot insert the plug into the outlet or if you are not sure whether the outlet is grounded, contact your local electrician.

The instrument should only be operated from a power source indicated on the product's electrical ratings label. If you have questions about the type of power source to use, contact your authorized Miltenyi Biotec service provider or local power company.

Make sure that the mains switch as well as the connector for the power cable are easily accessible and located as close to the operator of the instrument as possible. If it is necessary to disconnect the power supply, unplug the cable from the power outlet.

Servicing

Unless otherwise mentioned in the Miltenyi Biotec documentation, do not service the equipment yourself. Service needed on components inside the instrument must be carried out only by personnel trained and authorized by Miltenyi Biotec.

Chemical and biological hazards

The proper handling of the gentleMACS™ Tubes minimizes exposure to the samples. However, potential hazards to the operator still exist. The operator must choose and wear personal safety equipment in accordance with warnings and precautions for the used substances. Defective or inadequate safety equipment might endanger the operator.

- The gentleMACS Dissociator must be operated within a safety enclosure when hazardous or unknown materials are processed, to protect persons and the laboratory environment from contamination.

- All liquid and solid waste must be considered biohazardous and, therefore, must be handled taking universal laboratory precautions.
- All clinical samples must be considered potentially infectious.
- Toxic or corrosive chemicals, poisonous, acidic, or radioactive substances may be present.

- It is strictly prohibited for anyone other than the operator to handle contaminated accessories or contaminated parts of the gentleMACS Dissociator.

- After moving the gentleMACS Dissociator from a cold environment, for example, a cold room at 4 °C to room temperature, it is necessary to wait a minimum of 4 hours before setting the instrument into operation.
- Condensing liquid droplets can severely damage the instrument.

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The gentleMACS™ Dissociator is a benchtop instrument for the automated dissociation of tissues to obtain single-cell suspensions. Single-cell suspensions are a crucial prerequisite to achieve reliable results in cell separation, cell analysis, and cell culture. Moreover, the gentleMACS Dissociator can be used for the homogenization of cells and tissues to extract biomolecules for molecular biology or protein biochemistry applications.

Two types of gentleMACS Tubes have been specially designed: C Tubes for the gentle preparation of single-cell suspensions from tissues and M Tubes for the homogenization of tissues and cells. Both types of tubes are designed for safe and sterile sample handling.

A number of gentleMACS Programs have been specifically developed to offer standardized procedures for the gentle dissociation or homogenization of various tissues. A single sample or two samples in parallel can be processed.

A comprehensive gentleMACS Program overview and specific protocols for numerous applications are available for download at **www.gentlemacs.com**.

3 gentleMACS™ Products and Spare Parts

gentleMACS™ Products

Product	Description	Order no.
gentleMACS™ Starting Kit	gentleMACS Dissociator C Tubes (25 pieces) M Tubes (25 pieces) 1 Power cord 1 gentleMACS Dissociator user manual	130-093-235
C Tubes	For the dissociation of tissues to obtain single-cell suspensions. 25 tubes, single-packed	130-093-237
M Tubes	For the homogenization of tissues or cells to isolate subcellular material. 25 tubes, single-packed	130-093-236
	50 tubes per bag	130-093-458
M Tubes with Strainer	50 tubes per bag (strainer: 600 µm mesh)	130-094-392

For more information, visit www.miltenyibiotec.com

gentleMACS™ Spare Parts

Product	Description	Order no.
gentleMACS™ Sleeves	2 adapters for attachment of gentleMACS Tubes	130-093-464
gentleMACS™ Trough	Container for liquid waste	130-093-465
gentleMACS™ Safety Shield	Safety shield for protection	130-093-466
gentleMACS™ Program Card	Memory card containing gentleMACS Programs	130-093-462

Only use accessories recommended by Miltenyi Biotec. Failure to use recommended spare parts and disposables may compromise safe operation and may result in inaccurate results and reduced life time of the instrument. Miltenyi Biotec does not honor any warranty or accept any responsibility for damages resulting from the use of inappropriate accessories.

4 Storage and handling

The gentleMACS™ Dissociator must be stored and operated under dry conditions (non-condensing humidity). Storage temperature must be between $-15\text{ }^{\circ}\text{C}$ and $+70\text{ }^{\circ}\text{C}$. The instrument should be operated at a temperature between $+4\text{ }^{\circ}\text{C}$ and $+35\text{ }^{\circ}\text{C}$.

Please unpack the instrument carefully and check for completeness and for any damage that may have occurred during transportation. If there is any damage, contact your local Miltenyi Biotec representative or Miltenyi Biotec technical service. Please keep the original packaging material.

The gentleMACS™ Dissociator is designed exclusively for indoor use.

Components included in delivery

- 1 gentleMACS™ Dissociator
- C Tubes (25 pieces)
- M Tubes (25 pieces)
- 1 Power cord
- 1 gentleMACS Dissociator user manual

Installation

The gentleMACS™ Dissociator requires a flat surface with the dimensions 250 mm × 400 mm and a power connection. To operate the instrument, the keyboard area and the top of the instrument must be accessible. The gentleMACS Dissociator is a ready-to-use instrument.

6 Operating the gentleMACS™ Dissociator

General information

The gentleMACS™ Dissociator is used as a stand-alone, benchtop instrument. A number of pre-defined programs offer a variety of processes for tissue dissociation and tissue homogenization. The specially designed gentleMACS Tubes are required for processing the samples. The gentleMACS Dissociator must be used exclusively with these Tubes.


C Tubes are designed for the dissociation of tissues to obtain single-cell suspensions. M Tubes are designed for homogenization of tissues or cells to isolate subcellular material. A single sample or two samples in parallel can be processed.

For details on the sample preparation, see the respective gentleMACS Protocols, which are available at www.miltenyibiotec.com/protocols.

Below, please find brief instructions for the use of the instrument. For detailed information, see the next sections.

- Make sure that the gentleMACS Dissociator is switched ON. The mains switch is located at the back panel of the instrument.
- Pre-defined programs are provided by the internal gentleMACS Dissociator memory.
Alternatively: Insert the gentleMACS Program Card to run a program from the Program Card or to copy a program from the Program Card to the internal memory of the instrument.
- Transfer the samples into the Tubes. Make sure that the Tubes are securely closed.
- Make sure to properly install the Tubes vertically and upside down onto the Sleeves. Proper installation will ensure that the Tubes are held in position against rotational and axial forces.
- All functions of the gentleMACS Dissociator must be controlled via the display; see the following sections for details.
- Select one of the programs offered by the instrument memory or by the gentleMACS Program Card.
- Start the process by pressing the start button.
- The display indicates the program name and the remaining time (seconds) for the current program to finish.
- The end of the process is displayed for 5 seconds. Subsequently, the display indicates the program that has just been used. By pressing the start button, this program can be applied immediately to process the next sample.
- Remove the Tubes containing the sample by pulling it vertically from the Sleeve.

Layout of the keyboard

Panel buttons

Selection of one of the three modes that allow the user to select a program for sample processing, or copy, or delete programs

Selection of a program

Confirmation of a selection or function

Run or stop a process

gentleMACS™ Program Card

Insert the gentleMACS™ Program Card in the direction indicated by the arrow. The card will then automatically be selected to access programs.

If the gentleMACS Program Card is inserted, only the programs from the Program Card are available. The programs from the Program Card can be copied to the internal gentleMACS Dissociator memory. To select programs from the internal memory, remove the gentleMACS Program Card.

Function of the LED

The LED indicates the instrument status for convenient monitoring.

Green:	System ready
Green flashing:	Process running
Red:	System error
Red flashing:	Error during process

Display

Using the menu button, the operator can switch between select, copy, and delete modes.

Select mode

<Program name>
<Tube name>

Copy mode

Copy
<Program name>

Delete mode

Delete
<Program name>

Power ON

gentleMACS V1.xx
Miltenyi Biotec

Display during initialization.


<Program name>
<Tube name>

After initialization, the green LED is on and the gentleMACS™ Dissociator is ready for operation.

Select and run a program

To select a program from the internal gentleMACS™ Dissociator memory, press the menu button until the select menu is displayed. If the gentleMACS Program Card is inserted, only the programs from this card are available.


To select programs from the gentleMACS Dissociator memory, remove the gentleMACS Program Card. To copy programs from the gentleMACS Program Card to the gentleMACS Dissociator memory, see the section Copy a program.

6

Copy a program

To copy a program from the gentleMACS™ Program Card to the gentleMACS Dissociator memory, insert the gentleMACS Program Card. Press the menu button until the copy menu is displayed. If the gentleMACS Program Card is not inserted, a program from the gentleMACS Dissociator memory will be selected.

Select the desired program from the gentleMACS Program Card using the arrow buttons. If a program is copied within the internal memory, "Copy memory" is displayed instead of "Copy card".

Press **ok**.

Remove the gentleMACS Program Card.

Select the destination on the gentleMACS Dissociator memory using the arrow buttons. If there is already a program stored at the selected destination, the name of this program is displayed. If the destination is empty, the line remains blank.

Press **ok** to store the program or press **menu** to cancel and select a different program.


If there is already a program stored at the selected destination, confirm "overwrite?" by pressing **ok** or cancel by pressing **menu**.

The program is being copied.

If an error occurs, confirm the error message by pressing **stop**.

Delete a program

To delete a program from the gentleMACS™ Dissociator memory, press the menu button until the delete menu is displayed.

Emergency shutdown procedure

1. Press the **stop** button.
2. Press the mains switch at the back panel of the instrument.
3. Unplug the instrument.

7 Performance and limitations of use

- A maximum of 99 programs can be stored on the internal gentleMACS Dissociator memory.
- Editing programs on the gentleMACS Dissociator is not possible.
- If an error occurs, the gentleMACS Dissociator displays a message that needs to be confirmed by the operator. Press the **stop** button in order to continue operation.
- Both motors run at the same time, theoretically at identical speed and acceleration. However, the friction and resistance of the sample will decrease the effective speed. The gentleMACS™ Dissociator will display an error message if a pre-set tolerance has been exceeded during the process.
- If one of the motors is blocked due to high resistance of the sample in the tube, the instrument will try four times to get started again. The message "Error: Motor 1 / 2 blocked" appears and the red LED flashes. The operator must confirm the incident by pressing the **stop** button.

If the blockage can be resolved, the gentleMACS Dissociator continues the procedure from where it was interrupted. The message "Warning: x Retries" appears for 5 seconds during the incident and disappears automatically.

If the blockage cannot be resolved, the power of the motor is turned OFF.

Cleaning

In case of spillage or contamination, clean the instrument appropriately. For decontamination, use 70% ethanol or a 1% hypochlorite solution.

- Do not allow fluids or objects to enter the interior of the instrument.
- Make sure to clean the instrument exclusively with cleaning agents recommended by Miltenyi Biotec. Do not use acetone or any aggressive agents.
- Do not leave Tubes attached to the instrument while it is not in use.
- When processing infectious, radioactive, poisonous, or any other hazardous liquids, always abide by the necessary safety precautions.
- Before using organic solvents or aggressive chemicals in the dissociation or homogenization process, the compatibility with the gentleMACS™ Dissociator parts and Tubes must be tested.
- Large amounts of accidentally spilled liquid are automatically collected in the gentleMACS Trough. Use a pipette to aspirate the liquid before removing the Trough for cleaning and sterilizing.
- If servicing of the gentleMACS Dissociator is required, send the instrument back to a service provider authorized by Miltenyi Biotec. Make sure that all parts are decontaminated and sterilized. Use the original packaging material.


- Always use your personal safety equipment while removing or disassembling parts of the gentleMACS Dissociator for cleaning or sterilization.

Cleaning procedure

For details on disassembly and reassembly, see the picture below.

1. Remove the Safety Shield.
2. Remove the Sleeves from the instrument for cleaning and/or sterilization.
3. Remove the Trough for cleaning and/or sterilization.
4. The remaining surface of the gentleMACS™ Dissociator should be cleaned with a cloth soaked with a detergent-containing aqueous solution. Carefully wipe all parts.
5. Repeat the procedure with another cloth soaked with distilled water.
6. If the homogenized substances are not water-soluble, it is recommended to wipe all parts with ethanol or an alternative organic solvent. Contact Miltenyi Biotec technical support for further advice.
7. The removed parts should be washed in a laboratory dishwasher. It is recommended to use a cleaning agent containing detergent and alkaline cleaning solution.
8. If necessary, exchange the Sleeves and the Trough.
9. Reassemble the parts and make sure that the removed parts snap into position to ensure the correct positioning.

Disassembly for cleaning and reassembly

- Safety Shield, Sleeves, and Trough are removed by pulling precisely vertically with a slight force.
- The Trough can be detached easily by slightly pressing the rippled sides.
- Do not use any tools to detach the described parts.
- Reassemble the parts in reverse order.
- Correct assembly of all parts can be sensed by a slight snap as the part slips into its final position.


The gentleMACS™ Dissociator and the removable parts must be cleaned as described above. The lifetime of the Sleeves might be limited due to the mechanical strain during the homogenization procedure. Periodically, check the parts for wear.

Attention: Decontamination is within the responsibility of the operator of the instrument. Contaminated instruments or removable parts must not be given to any person. After detachment and reassembly of any removable part, it is necessary to check that the part still functions properly.

Applied materials

- 1 Housing middle part: PA 66/6 (with 30% glass fibers)
- 2 Housing front part: PA 66/6 (with 30% glass fibers)
- 3 Housing rear part: PA 66/6 (with 30% glass fibers)
- 4 Keyboard: Polyester film
- 5 gentleMACS™ Safety Shield: PP
- 6 gentleMACS™ Trough: PA 66/6 (with 30% glass fibers)
- 7 gentleMACS™ Sleeve: PA 66/6 (with 30% glass fibers)

Internal quality control

During the lifetime of the instrument no calibration is necessary. To check the quality of the instrument, the gentleMACS™ Dissociator must be sent to Miltenyi Biotec.

Disposal information

The gentleMACS™ Dissociator shall be disposed of according to the biosafety level of its use. Potentially contaminated parts must be sterilized before being disposed of, according to local laboratory regulations.

Attention: Decontamination is within the responsibility of the operator of the instrument. Contaminated instruments or contaminated parts of the instrument must not be sent or given to any person.

Error	Source of error	User action
Liquid droplets in Trough	Tube is not closed properly.	Close Tube tightly beyond the first resistance.
	Content of Tube is jammed between cap and Tube.	The top area of cap and Tube must be perfectly clean for proper functionality of the seal.
	Tube seal is leaky.	Replace Tube. Clean the instrument as described in chapter 8.
Tube vibrates during run	Sleeve is not installed correctly.	Ensure correct installation of the Sleeve. Make sure that the Sleeve snaps into position.
	Inappropriate sample material.	Refer to the respective datasheet for technical specifications of the gentleMACS™ Tubes.
Tube is not tightly held in position during run	Tube does not fit properly into the snap-fit of the Trough.	Ensure correct installation of the Tube.
	Tube is not closed properly.	Close Tube tightly beyond the first resistance.
gentleMACS™ Dissociator is without power	Power cord is not connected.	Check power cord.
	Instrument is not switched ON.	Check mains switch at the rear panel of the instrument.
	Fuse is blown.	Check fuses at the rear panel of the instrument. Replace fuses only with those of the same type and rating.
gentleMACS Dissociator does not read the Program Card	Dust deposit on Program Card, or Program Card is damaged.	Check Program Card.

If any of the problems persist, please contact our technical support team (see chapter 12).

10 Technical specifications

Technical specifications		
Size	Length × Width × Height	347 mm × 215 mm × 195 mm
Weight		4.7 kg
Main power	Power consumption Supply voltage Fuse	max. 320 VA 100–240 VAC, 50–60 Hz 2 × T4AH250V
Motor	Rotational speed Rotation Rotational acceleration Torque	200–4,000 rpm clockwise or counter clockwise 100 r/s ² minimum 150 mNm. Overload restriction over the range of 200–4,000 rpm
Conditions of operation	Working temperature: Humidity Altitude Storage temperature	+4 °C to +35 °C 15–80%, relative, non-condensing maximum 2000 m –15 °C to +70 °C
USB interface	For service purposes only	Pin 1: USB +5 V Pin 2: USB – Pin 3: USB + Pin 4: USB GND

The gentleMACS Dissociator is labeled as a protection class I device and must be plugged into a grounded power outlet.

Supply voltage fluctuations up to +/- 10% of the nominal voltage. Transient overvoltages present on the mains supply: category II. The instrument is suitable for rated pollution degree 2. Emission sound pressure level at workstation < 70 dB(A). Design and specifications are subject to change without notice.

The gentleMACS™ Dissociator is designed to meet the requirements of the following standard:

- FCC 47 CFR Ch.1 part 15, class B (electromagnetic compatibility)

The gentleMACS™ Dissociator has been investigated by Underwriters Laboratories in accordance with the following standards:

- UL 61010-1
- CAN/CSA-C22.2 No. 61010-1
- IEC 61010-2-051.

For other safety considerations, refer to the product label, or visit www.miltenyibiotec.com.

Except as stated in a specific warranty statement, which may accompany this product, or unless otherwise agreed in writing by an authorized representative of Miltenyi Biotec, Miltenyi Biotec's warranty, if any, with respect to this product is subject to the terms and conditions of sale of the company within the Miltenyi Biotec group which supplied the product. Copies are available on request or at www.miltenyibiotec.com. Terms and conditions of sale may vary by country and region. Nothing herein should be construed as constituting an additional warranty.

Miltenyi Biotec's warranty for this product only covers product issues caused by defects in material or workmanship during ordinary use; it does not cover product issues caused by any other reason, including but not limited to product issues due to use in a manner other than specifically described in this manual, for example: inappropriate or improper use; incorrect assembly or installation by an operator or a third party; reasonable wear-and-tear; negligent or incorrect handling, servicing, or maintenance; non-adherence to the operating instructions; unauthorized modification of or to any part of this product; or use of inappropriate accessories or work materials. Miltenyi Biotec's warranty does not cover products sold AS IS or WITH ALL FAULTS or consumables. Nothing herein should be construed as constituting an additional warranty.

Miltenyi Biotec must be informed immediately, if a claim is made under such warranty. If a material or manufacturing defect occurs within the warranty period, Miltenyi Biotec will take the appropriate steps to restore the full usability of the unit.

Limitation on damages:

Miltenyi Biotec shall not be liable for any incidental or consequential damages for breach of any express or implied warranty or condition on this product.

Some states or jurisdictions do not allow the exclusion or limitation of incidental or consequential damages, so the above limitations or exclusions may not apply to you. This warranty statement gives you specific legal rights and you may have other rights, which vary from state to state or jurisdiction to jurisdiction.

12 Technical service

For technical questions, please contact your local distributor or our technical support team:

Miltenyi Biotec GmbH
Friedrich-Ebert-Str. 68
51429 Bergisch Gladbach
Germany
Phone +49 2204 8306-830
Fax +49 2204 8306-89
macstec@miltenyibiotec.de
www.miltenyibiotec.com

Name and address of the manufacturer:

Miltenyi Biotec GmbH
Friedrich-Ebert-Str. 68
51429 Bergisch Gladbach
Germany

The declaration of conformity refers to the machinery identified as follows:

Description: laboratory equipment for mixing and stirring
Function: dissociation of tissues to obtain single-cell suspensions or homogenization of tissues
Model: gentleMACS™ Dissociator

We hereby declare that the machinery complies with all essential requirements of the following directives:

2004/108/EC (electromagnetic compatibility)
2006/42/EC (machinery)

and is in conformity with the following harmonized standards:

EN 61010-1:2001
EN 61010-2-051:2003
EN 61326-1:2006

Person authorized to compile the relevant technical documentation:

Stefan Miltenyi
Miltenyi Biotec GmbH
Friedrich-Ebert-Str. 68
51429 Bergisch Gladbach
Germany

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Miltenyi Biotec provides products and services worldwide. Visit www.miltenyibiotec.com/local to find your nearest Miltenyi Biotec contact.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. MACS is a registered trademark and gentleMACS is a trademark of Miltenyi Biotec GmbH. Copyright © 2010 Miltenyi Biotec GmbH. All rights reserved.

Miltenyi Biotec

gentleMACS™ Dissociator User manual

Addendum

Important information

Please adhere to the following precaution:
The instrument must be used only with properly installed gentleMACS™ Tubes whose lids are tightly closed. Do not operate the instrument with gentleMACS Tube lids only.

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Miltenyi Biotec provides products and services worldwide. Visit www.miltenyibiotec.com/local to find your nearest Miltenyi Biotec contact.

MACS is a registered trademark and gentleMACS is a trademark of Miltenyi Biotec GmbH.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use.
Copyright © 2010 Miltenyi Biotec GmbH. All rights reserved.